



Everett
1300 44th Street SE
Everett, WA 98203
Phone 425-339-2559
Fax 425-339-1583

Bellingham
3106 NW Ave.
Bellingham, WA
98225
Phone 360-734-0298
Fax 360-734-9679

Mt. Vernon
1600 Continental
PL, STE 105
Mt. Vernon, WA
98273
Phone 360-416-6505
Fax 360-416-8241

Edmonds
7631 212th Street
SW, STE 102A
Edmonds, WA
98026
Phone 425-776-1247
Fax 425-776-8045

Sequim
566 N Fifth Ave
Sequim, WA 98382
Phone 360-797-1001
Fax -360-797-1003

Spokane
101 E Hastings Rd,
STE J
Spokane, WA 99218
Phone 509-620-9362
Fax 509-505-0912

**INSTRUCTIONS FOR
CARE OF PROSTHETIC SOCKS**

Prosthetic socks are used by amputees to maintain proper positioning of the residual limb inside the prosthesis. Additionally, the prosthetic sock protects and cushions the skin, absorbs perspiration, and aids in reducing friction between the prosthesis and the residual limb.

**PROSTHETIC SOCKS SHOULD BE CHANGED DAILY AND WASHED IN WARM OR COOL WATER
WITH MILD SOAP**

We advise that the socks be washed separately from other laundry, as dyes from other clothes can become imbedded in the prosthetic sock and may irritate the residual limb.

**PROSTHETIC SOCKS ARE WOVEN ESPECIALLY FOR THEIR INTENDED USE AND ARE AVAILABLE IN
VARIOUS SIZES AND THICKNESSES**

Your prosthetist will determine the size that is correct and generally will supply you with:

- 1. A nylon sheath (worn next to the skin)
- 2. A 1-ply poly-soft sock (stretches, white sock)
- 3. A 3-ply wool sock
- 4. A 5-ply wool sock

With your first prosthesis, you will be supplied with approximately six prosthetic socks. If you need more socks, please contact your prosthetist.

The number of ply relates to the thickness of the prosthetic sock and aids in maintaining proper fit of the prosthesis. A 5-ply is made out of five threads of yarn, whereas, a 3-ply sock is made out of 3 threads of yarn. When these two socks are used together, there are a total of eight (8) plies of sock on the residual limb. It is important to remember that as the prosthetic sock becomes worn, the thickness of the sock changes, i.e. an old 5-ply sock would perhaps equal a new 3-ply sock.

It is also important to remember that the proper amount of socks that you should be wearing is what feels comfortable to you. If the prosthesis feels too loose, try adding a ply of socks. If it feels too tight, try decreasing the number of ply. If the prosthesis still does not feel comfortable to you, contact your prosthetist as additional padding in the socket may be necessary.

Prosthetic socks are available in a wide variety of thicknesses, lengths, and materials. If you have any questions on the selection and care of prosthetic socks, please contact your prosthetist.

If you have any problems with the use of your device; discontinue using it immediately and contact Cornerstone.

Signature _____

Date _____